Opis produktu

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Wydanie: [10.2021]

Model: DRQM6J

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Edition: [10.2021]

Model: DRQM6J

- Advanced wrinkle prevention
- Better skin elasticity and firmness
- Spa effect at home

Product specification:

Dimensions: 15,2x5,8x1,6 cm

Weight: 77 g

Vibration frequency: 6000/min

Composition: aluminum alloy, zinc alloy, ABS

Power supply: 1x AA battery (not included)

Colors: silver, pink

Operating temperature: 0-40°C

Product complies with standards:

EN 55014-1: 2017+A11: 2020

EN 55014-2: 2015

Qmed 3D Massager has been developed with the aim to combine best practices in the traditional treatment using touch therapy, modern engineering and the art of creating beautiful items. The history of face massage started in antiquity, when cosmetology – art of beautifying – maintained a prominent position in medicine, usually accessible only for ruling royal families or noble wealthiest citizens. Queen Cleopatra – famous for her beauty – was everyday massaged and anointed with oils. In Japanese tradition beauty is a sign of health condition, interpreted as harmony between physical, mental and spiritual zones. The aim of the massage is to support this harmony and to maintain natural beauty. Diamond-grade cut of the heads and their arrangement in V-shape at the angle of 70°, perfectly fits the product to body's line and allows for massage comparable to spa treatment. 3D Face massager is recommended for contouring face, sharpening the cheekbones, reducing swelling and dark circles under the eyes, improving the oval of the face, firming the skin of the neck and décolletage. Positive effects are visible with a 10 minutes daily massage applying favorite oil or serum, as well as using the device on the sheet mask increasing absorption of beneficial ingredients. Waterproof

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construction of the product allows for massage during relaxing bath.

Important! Although metal does not have properties favoring bacterial growth, remember to massage with clean device and disinfect it after each use.

Do not apply too much pressure, gentle movements are sufficient. Pressing down the massager to the body does not increase the results. Do the massage from the center of the face outwards, repeating each movement several times. Do not use the product if you have skin problems, discontinuity of epidermis or allergy to metal.

Example of massage direction:

The massager may be used for professional purposes or for personal home use. The device does not generate an awkward noise, which is particularly important if relaxation is the aim.

ATTENTION: Always follow the basic precautions according to the information included below and regarding safety measures when using the device. Read the entire instructions carefully before use. Proper use of this device is recommended for best results. Improper use of this device could result in injury or health complications. Keep this user manual for future reference. The manufacturer reserves the right to alter the specifications of the product without prior notification.

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Warnings:

1. Keep safe from falling down. An impact may damage the product.
2. Do not use the device in a manner inconsistent with the instructions.
3. Carefully check your device before each use.
4. Do not make any private modifications to the equipment.
5. After using and cleaning the massager, wipe the surface with a soft cloth. Make sure the product is dry before putting it back to its packaging. Do not use solvents, toilet cleaners, sharp brushes, hard objects and caustic cleaners. Do not attempt to clean internal surfaces. Disassembling the device to clean internal mechanisms could result in device failure.
6. Storing the massager in areas with very high humidity (baths, laundry, etc.) or in the open air as well as lack of hygiene causes rapid loss of operating parameters and aesthetic of the product.
7. It is unacceptable to use the massager after discovering any visible or detectable damages that may cause an accident.
8. Do not let children or persons who are not capable of handling equipment operate the product.
9. Do not use the product for any purpose, other than massage. This may result in an accident or injury.
10. Please keep the product out of reach of children.

Maintenance:

1. Check if there is any damage or wear by visual inspection.
2. Be sure to remove the AC plug from the electrical outlet before cleaning the massager.
3. Wipe the surface of the device with a slightly damp towel and dry with a soft cloth. Make sure the massager is dry before putting them back in their packaging. Do not use solvents, toilet cleaners, sharp brushes and hard objects and caustic cleaners. Do not attempt to clean internal surfaces. Disassembling the device to clean internal mechanisms could result in device failure.
4. When storing or traveling for a long time, separate the device from the battery and place it in the device's original case.

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5. It is extremely important to keep all the massager components clean to maintain the product parameters.

6. Storing the massager in areas with very high humidity (baths, laundry, etc.) or in the open air as well as lack of hygiene causes rapid loss of operating parameters and aesthetic of the product.

7. It is unacceptable to use the massager after discovering any visible or detectable damages or when the wear of some parts may cause an accident. All the repairs need to be performed by skilled personnel and during the warranty period only by the manufacturer's service.

8. Temperature for storage of the massager shall be higher than 0°C and relative moisture ≤60%.

Disclaimer:

Please read the instruction manual carefully before using this equipment. Use this equipment correctly in accordance with the instructions. MDH Sp. z o. o. will not be liable for any damage or loss to the customer caused by the incorrect use of this equipment as explained.

The product complies with the requirements of the standards for electrical and electronic devices. It should not be disposed of with other household waste. The user is responsible for handing over the device at the end of its service life to the appropriate collection points, under penalty of sanctions provided for by the regulations in force in the matter of waste disposal. For more information on the collection systems available, please contact your local waste disposal center.

The package contains 1 AA battery-powered product, the battery is not included.

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